

Product Datasheet

Human Pregnancy-specific Beta-1-glycoprotein 9 ELISA Kit (orb1209546)

Description	Human Pregnancy-specific Beta-1-glycoprotein 9, PSG9 ELISA Kit
Reactivity	Human
Range	0.5-300ug/mL
Tested Applications	ELISA
Target	PSG9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	48,272
Uniprot ID	Q00887
Solubility (25°C)	320µg/mL
Sensitivity	0.25µg/mL
Entrez	5678
Expiration Date	Please enquire.