

Product Datasheet

Human Platelet Factor 4 Variant ELISA Kit (orb1209556)

Description	Human Platelet Factor 4 Variant, PF4V1 ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	PF4V1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,553
Uniprot ID	P10720
Solubility (25°C)	96ng/mL
Sensitivity	0.143ng/mL
Entrez	5197
Expiration Date	Please enquire.