

Product Datasheet

Human Pikachurin ELISA Kit (orb1209559)

Description	Human Pikachurin, EGFLAM ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	EGFLAM
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	111,271
Uniprot ID	Q63HQ2
Solubility (25°C)	4800ng/L
Sensitivity	9.26ng/L
Entrez	133584
Expiration Date	Please enquire.