

Product Datasheet

Human Phosphatidylethanolamine-binding Protein 1 ELISA Kit (orb1209570)

Description	Human Phosphatidylethanolamine-binding Protein 1, PEBP1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	PEBP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,057
Uniprot ID	P30086
Solubility (25°C)	24ng/mL
Sensitivity	0.029ng/mL
Entrez	5037
Expiration Date	Please enquire.