

Product Datasheet

Human Peptidase Inhibitor 16 ELISA Kit (orb1209579)

Description	Human Peptidase Inhibitor 16, PI16 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	PI16
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	49,471
Uniprot ID	Q6UXB8
Solubility (25°C)	24ng/mL
Sensitivity	0.026ng/mL
Entrez	221476
Expiration Date	Please enquire.