

Product Datasheet

Human Oligodendrocyte-Myelin Glycoprotein ELISA Kit (orb1209586)

Description	Human Oligodendrocyte-Myelin Glycoprotein, OMG ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	OMG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	49,608
Uniprot ID	P23515
Solubility (25°C)	24ng/mL
Sensitivity	0.029ng/mL
Entrez	4974
Expiration Date	Please enquire.