

Product Datasheet

Human Odontogenic Ameloblast-associated Protein ELISA Kit (orb1209589)

Description	Human Odontogenic Ameloblast-associated Protein, ODAM ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	ODAM
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,777
Uniprot ID	A1E959
Solubility (25°C)	24ng/mL
Sensitivity	0.024ng/mL
Entrez	54959
Expiration Date	Please enquire.