

Product Datasheet

Human Nuclear Receptor Subfamily 4 Group A Member 2 ELISA Kit (orb1209593)

Description	Human Nuclear Receptor Subfamily 4 Group A Member 2, NR4A2 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	NR4A2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	66,591
Uniprot ID	P43354
Solubility (25°C)	4800ng/L
Sensitivity	9.19ng/L
Entrez	4929
Expiration Date	Please enquire.