

Product Datasheet

Human Neuroendocrine Protein 7B2 ELISA Kit (orb1209607)

Description	Human Neuroendocrine Protein 7B2, SCG5 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	SCG5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,730
Uniprot ID	P05408
Solubility (25°C)	24ng/mL
Sensitivity	0.027ng/mL
Entrez	6447
Expiration Date	Please enquire.