

Product Datasheet

Human Musculoskeletal Embryonic Nuclear Protein 1 ELISA Kit (orb1209620)

Description	Human Musculoskeletal Embryonic Nuclear Protein 1, MUSTN1 ELISA Kit
Reactivity	Human
Range	20-3500ng/L
Tested Applications	ELISA
Target	MUSTN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	8,911
Uniprot ID	Q8IVN3
Solubility (25°C)	4000ng/L
Sensitivity	8.98ng/L
Entrez	389125
Expiration Date	Please enquire.