

Product Datasheet

Human Melanoma-associated Antigen D2 ELISA Kit (orb1209638)

Description	Human Melanoma-associated Antigen D2, MAGED2 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	MAGED2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	64,954
Uniprot ID	Q9UNF1
Solubility (25°C)	4800ng/L
Sensitivity	9.82ng/L
Entrez	10916
Expiration Date	Please enquire.