

Product Datasheet

Human Matrix Metalloproteinase-28 ELISA Kit (orb1209640)

Description	Human Matrix Metalloproteinase-28, MMP28 ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	MMP28
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	58,939
Uniprot ID	Q9H239
Solubility (25°C)	96ng/mL
Sensitivity	0.156ng/mL
Entrez	79148
Expiration Date	Please enquire.