

Product Datasheet

Human Lysosomal Protective Protein ELISA Kit (orb1209655)

Description	Human Lysosomal Protective Protein, CTSA ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	CTSA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,466
Uniprot ID	P10619
Solubility (25°C)	2400ng/L
Sensitivity	4.95ng/L
Entrez	5476
Expiration Date	Please enquire.