

Product Datasheet

Human Lysosomal Alpha-glucosidase ELISA Kit (orb1209657)

Description	Human Lysosomal Alpha-glucosidase, GAA ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	GAA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	105,324
Uniprot ID	P10253
Solubility (25°C)	3200ng/L
Sensitivity	7.12ng/L
Entrez	2548
Expiration Date	Please enquire.