

Product Datasheet

Human Lysosomal Acid Lipase Cholesteryl Ester Hydrolase ELISA Kit (orb1209658)

Description	Human Lysosomal Acid Lipase Cholesteryl Ester Hydrolase, LIPA ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	LIPA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	45,419
Uniprot ID	P38571
Solubility (25°C)	9600ng/L
Sensitivity	14.22ng/L
Entrez	3988
Expiration Date	Please enquire.