

Product Datasheet

Human Lysine-specific Demethylase 5A ELISA Kit (orb1209659)

Description	Human Lysine-specific Demethylase 5A, KDM5A ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	KDM5A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	192,095
Uniprot ID	P29375
Solubility (25°C)	3200ng/L
Sensitivity	7.03ng/L
Entrez	5927
Expiration Date	Please enquire.