

Product Datasheet

Human Lymphocyte Antigen 96 ELISA Kit (orb1209661)

Description	Human Lymphocyte Antigen 96, LY96 ELISA Kit
Reactivity	Human
Range	50-10000ng/L
Tested Applications	ELISA
Target	LY96
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,546
Uniprot ID	Q9Y6Y9
Solubility (25°C)	16000ng/L
Sensitivity	24.36ng/L
Entrez	23643
Expiration Date	Please enquire.