

Product Datasheet

Human Ly6 Plaur Domain-containing Protein 3 ELISA Kit (orb1209665)

Description	Human Ly6 Plaur Domain-containing Protein 3, LYPD3 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	LYPD3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,971
Uniprot ID	O95274
Solubility (25°C)	4800ng/L
Sensitivity	9.11ng/L
Entrez	27076
Expiration Date	Please enquire.