

## Product Datasheet

# Human Ly6,Plaur Domain-containing Protein 1 ELISA Kit (orb1209666)

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Description</b>         | Human Ly6,Plaur Domain-containing Protein 1, LYPD1 ELISA Kit |
| <b>Reactivity</b>          | Human                                                        |
| <b>Range</b>               | 15-3000ng/L                                                  |
| <b>Tested Applications</b> | ELISA                                                        |
| <b>Target</b>              | LYPD1                                                        |
| <b>Storage</b>             | 2-8°C                                                        |
| <b>Note</b>                | For research use only                                        |
| <b>Application notes</b>   | Sandwich                                                     |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                       |
| <b>MW</b>                  | 15,240                                                       |
| <b>Uniprot ID</b>          | <b>Q8N2G4</b>                                                |
| <b>Solubility (25°C)</b>   | 3200ng/L                                                     |
| <b>Sensitivity</b>         | 7.11ng/L                                                     |
| <b>Entrez</b>              | <b>116372</b>                                                |
| <b>Expiration Date</b>     | Please enquire.                                              |