

Product Datasheet

Human Lrp2-binding Protein ELISA Kit (orb1209668)

Description	Human Lrp2-binding Protein, LRP2BP ELISA Kit
Reactivity	Human
Range	50-100000ng/L
Tested Applications	ELISA
Target	LRP2BP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,780
Uniprot ID	Q9P2M1
Solubility (25°C)	12800ng/L
Sensitivity	23.44ng/L
Entrez	55805
Expiration Date	Please enquire.