

Product Datasheet

Human Low Density Lipoprotein Receptor Adapter Protein 1 ELISA Kit (orb1209669)

Description	Human Low Density Lipoprotein Receptor Adapter Protein 1, LDLRAP1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	LDLRAP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,885
Uniprot ID	Q5SW96
Solubility (25°C)	2400ng/L
Sensitivity	4.88ng/L
Entrez	26119
Expiration Date	Please enquire.