

Product Datasheet

Human Loricrin ELISA Kit (orb1209671)

Description	Human Loricrin, LOR ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	LOR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,761
Uniprot ID	P23490
Solubility (25°C)	3200ng/L
Sensitivity	7.86ng/L
Entrez	4014
Expiration Date	Please enquire.