

Product Datasheet

Human L-lactate Dehydrogenase B Chain ELISA Kit (orb1209674)

Description	Human L-lactate Dehydrogenase B Chain, LDHB ELISA Kit
Reactivity	Human
Range	15-3000mU/mL
Tested Applications	ELISA
Target	LDHB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,638
Uniprot ID	P07195
Solubility (25°C)	3200mU/mL
Sensitivity	7.22mU/mL
Entrez	3945
Expiration Date	Please enquire.