

Product Datasheet

Human Liver-expressed Antimicrobial Peptide 2 ELISA Kit (orb1209675)

Description	Human Liver-expressed Antimicrobial Peptide 2, LEAP2 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	LEAP2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	8,814
Uniprot ID	Q969E1
Solubility (25°C)	48ng/ml
Sensitivity	0.058ng/ml
Entrez	116842
Expiration Date	Please enquire.