

Product Datasheet

Human Lipopolysaccharide-Responsive and Beige-like Anchor Protein ELISA Kit (orb1209677)

Description	Human Lipopolysaccharide-Responsive and Beige-like Anchor Protein, LRBA ELISA Kit
Reactivity	Human
Range	70-15000ng/L
Tested Applications	ELISA
Target	LRBA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	319,108
Uniprot ID	P50851
Solubility (25°C)	16000ng/L
Sensitivity	34.25ng/L
Entrez	987
Expiration Date	Please enquire.