

Product Datasheet

Human Leucine-rich Repeat Transmembrane Protein Flrt3 ELISA Kit (orb1209681)

Description	Human Leucine-rich Repeat Transmembrane Protein Flrt3, FLRT3 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	FLRT3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	73,004
Uniprot ID	Q9NZU0
Solubility (25°C)	6400ng/L
Sensitivity	9.72ng/L
Entrez	23767
Expiration Date	Please enquire.