

Product Datasheet

Human Leucine-rich Repeat Lgi Family Member 3 ELISA Kit (orb1209682)

Description	Human Leucine-rich Repeat Lgi Family Member 3, LGI3 ELISA Kit
Reactivity	Human
Range	0.1-12ng/mL
Tested Applications	ELISA
Target	LGI3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	61,704
Uniprot ID	Q8N145
Solubility (25°C)	24ng/mL
Sensitivity	0.022ng/mL
Entrez	203190
Expiration Date	Please enquire.