

Product Datasheet

Human IgG Receptor Fc γ Rn Large Subunit P51 ELISA Kit (orb1209717)

Description	Human IgG Receptor Fc γ Rn Large Subunit P51, FCGR1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	FCGR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,743
Uniprot ID	P55899
Solubility (25°C)	2560ng/L
Sensitivity	5.09ng/L
Entrez	2217
Expiration Date	Please enquire.