

Product Datasheet

Human Heat Shock Protein 75 kDa ELISA Kit (orb1209740)

Description	Human Heat Shock Protein 75 kDa, Mitochondrial, HSP75 ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	TRAP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	80,110
Uniprot ID	Q12931
Solubility (25°C)	96ng/mL
Sensitivity	0.14ng/mL
Entrez	10131
Expiration Date	Please enquire.