

Product Datasheet

Human Growth Differentiation Factor 6 ELISA Kit (orb1209746)

Description	Human Growth Differentiation Factor 6, GDF6 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	GDF6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,662
Uniprot ID	Q6KF10
Solubility (25°C)	4800ng/L
Sensitivity	9.54ng/L
Entrez	392255
Expiration Date	Please enquire.