

Product Datasheet

Human Granulocyte-macrophage Colony-stimulating Factor Receptor Subunit Alpha ELISA Kit (orb1209748)

Description	Human Granulocyte-macrophage Colony-stimulating Factor Receptor Subunit Alpha, CSF2RA ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	CSF2RA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,207
Uniprot ID	P15509
Solubility (25°C)	2400ng/L
Sensitivity	5.37ng/L
Entrez	1438
Expiration Date	Please enquire.