

Product Datasheet

Human Gpi-linked Nad ELISA Kit (orb1209749)

Description	Human Gpi-linked Nad, ART1 ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	ART1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,335
Uniprot ID	P52961
Solubility (25°C)	96ng/mL
Sensitivity	0.14ng/mL
Entrez	417
Expiration Date	Please enquire.