

Product Datasheet

Human Glia-derived Nexin ELISA Kit (orb1209758)

Description	Human Glia-derived Nexin, SERPINE2 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	SERPINE2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,002
Uniprot ID	P07093
Solubility (25°C)	48ng/ml
Sensitivity	0.052ng/ml
Entrez	5270
Expiration Date	Please enquire.