

Product Datasheet

Human G1 S-specific Cyclin-E1 ELISA Kit (orb1209760)

Description	Human G1 S-specific Cyclin-E1, CCNE1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CCNE1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,077
Uniprot ID	P24864
Solubility (25°C)	3200ng/L
Sensitivity	7.56ng/L
Entrez	898
Expiration Date	Please enquire.