

Product Datasheet

Human Fatty Acid-binding Protein 12 ELISA Kit (orb1209769)

Description	Human Fatty Acid-binding Protein 12, FABP12 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	FABP12
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,565
Uniprot ID	A6NFH5
Solubility (25°C)	25.6ng/mL
Sensitivity	0.023ng/mL
Entrez	646486
Expiration Date	Please enquire.