

## Product Datasheet

### Human Elav-like Protein 1 ELISA Kit (orb1209785)

|                            |                                                |
|----------------------------|------------------------------------------------|
| <b>Description</b>         | Human Elav-like Protein 1,<br>ELAVL1 ELISA Kit |
| <b>Reactivity</b>          | Human                                          |
| <b>Range</b>               | 15-3000ng/L                                    |
| <b>Tested Applications</b> | ELISA                                          |
| <b>Target</b>              | ELAVL1                                         |
| <b>Storage</b>             | 2-8°C                                          |
| <b>Note</b>                | For research use only                          |
| <b>Application notes</b>   | Sandwich                                       |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates         |
| <b>MW</b>                  | 36,092                                         |
| <b>Uniprot ID</b>          | <b>Q15717</b>                                  |
| <b>Solubility (25°C)</b>   | 3200ng/L                                       |
| <b>Sensitivity</b>         | 7.23ng/L                                       |
| <b>Entrez</b>              | <b>1994</b>                                    |
| <b>Expiration Date</b>     | Please enquire.                                |