

Product Datasheet

Human EF-Hand Domain-containing Protein D2 ELISA Kit (orb1209787)

Description	Human EF-Hand Domain-containing Protein D2, EFHD2 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	EFHD2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,697
Uniprot ID	Q96C19
Solubility (25°C)	2400ng/L
Sensitivity	5.49ng/L
Entrez	79180
Expiration Date	Please enquire.