

Product Datasheet

Human E3 Ubiquitin-protein Ligase Trim63 ELISA Kit (orb1209795)

Description	Human E3 Ubiquitin-protein Ligase Trim63, TRIM63 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	TRIM63
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,248
Uniprot ID	Q969Q1
Solubility (25°C)	1600ng/L
Sensitivity	3.13ng/L
Entrez	84676
Expiration Date	Please enquire.