

## Product Datasheet

### Human Dualoxidase 1 ELISA Kit (orb1209804)

|                            |                                         |
|----------------------------|-----------------------------------------|
| <b>Description</b>         | Human Dualoxidase 1,<br>DUOX1 ELISA Kit |
| <b>Reactivity</b>          | Human                                   |
| <b>Range</b>               | 15-3000ng/L                             |
| <b>Tested Applications</b> | ELISA                                   |
| <b>Target</b>              | DUOX1                                   |
| <b>Storage</b>             | 2-8°C                                   |
| <b>Note</b>                | For research use only                   |
| <b>Application notes</b>   | Sandwich                                |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates  |
| <b>MW</b>                  | 177,235                                 |
| <b>Uniprot ID</b>          | <b>Q9NRD9</b>                           |
| <b>Solubility (25°C)</b>   | 3200ng/L                                |
| <b>Sensitivity</b>         | 6.87ng/L                                |
| <b>Entrez</b>              | <b>53905</b>                            |
| <b>Expiration Date</b>     | Please enquire.                         |