

Product Datasheet

Human DNA Damage-inducible Transcript 4 Protein ELISA Kit (orb1209806)

Description	Human DNA Damage-inducible Transcript 4 Protein, DDIT4 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	DDIT4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,371
Uniprot ID	Q9NX09
Solubility (25°C)	24ng/ml
Sensitivity	0.018ng/ml
Entrez	54541
Expiration Date	Please enquire.