

Product Datasheet

Human Cytokine-like Protein 1 ELISA Kit (orb1209810)

Description	Human Cytokine-like Protein 1, CYTL1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	CYTL1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,577
Uniprot ID	Q9NRR1
Solubility (25°C)	24ng/ml
Sensitivity	0.021ng/ml
Entrez	54360
Expiration Date	Please enquire.