

Product Datasheet

Human Cysteinyl Leukotriene Receptor 1 ELISA Kit (orb1209815)

Description	Human Cysteinyl Leukotriene Receptor 1, CYSLTR1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CYSLTR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,541
Uniprot ID	Q9Y271
Solubility (25°C)	3200ng/L
Sensitivity	7.16ng/L
Entrez	10800
Expiration Date	Please enquire.