

Product Datasheet

Human Cst Complex Subunit Stn1 ELISA Kit (orb1209827)

Description	Human Cst Complex Subunit Stn1, OBFC1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	OBFC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,119
Uniprot ID	Q9H668
Solubility (25°C)	2400ng/L
Sensitivity	5.62ng/L
Entrez	79991
Expiration Date	Please enquire.