

Product Datasheet

Human Cd48 Antigen ELISA Kit (orb1209847)

Description	Human Cd48 Antigen, CD48 ELISA Kit
Reactivity	Human
Range	3-900ng/mL
Tested Applications	ELISA
Target	CD48
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,683
Uniprot ID	P09326
Solubility (25°C)	960ng/mL
Sensitivity	1.37ng/mL
Entrez	962
Expiration Date	Please enquire.