

Product Datasheet

Human Cd2-associated Protein ELISA Kit (orb1209848)

Description	Human Cd2-associated Protein, CD2AP ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	CD2AP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	71,451
Uniprot ID	Q9Y5K6
Solubility (25°C)	2400ng/L
Sensitivity	5.49ng/L
Entrez	23607
Expiration Date	Please enquire.