

Product Datasheet

Human Cd109 Antigen ELISA Kit (orb1209853)

Description	Human Cd109 Antigen, CD109 ELISA Kit
Reactivity	Human
Range	50-10000ng/L
Tested Applications	ELISA
Target	CD109
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	161,689
Uniprot ID	Q6YHK3
Solubility (25°C)	12800ng/L
Sensitivity	27.64ng/L
Entrez	135228
Expiration Date	Please enquire.