

Product Datasheet

Human C-C Chemokine Receptor Type 6 ELISA Kit (orb1209856)

Description	Human C-C Chemokine Receptor Type 6, CCR6 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	CCR6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,494
Uniprot ID	P51684
Solubility (25°C)	16ng/ml
Sensitivity	0.03ng/ml
Entrez	1235
Expiration Date	Please enquire.