

Product Datasheet

Human Catestatin ELISA Kit (orb1209857)

Description	Human Catestatin ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	catestatin
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	12,736
Uniprot ID	P61278
Solubility (25°C)	48ng/ml
Sensitivity	0.046ng/ml
Entrez	6750
Expiration Date	Please enquire.