

Product Datasheet

Human Adam Dec1 ELISA Kit (orb1209879)

Description	Human Adam Dec1, ADAMDEC1 ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	ADAMDEC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,775
Uniprot ID	015204
Solubility (25°C)	9600ng/L
Sensitivity	14.94ng/L
Entrez	27299
Expiration Date	Please enquire.