

Product Datasheet

Human Human 6 phosphofructo 2 kinase/fructose 2 ELISA Kit (orb1209882)

Description	Human Human 6 phosphofructo 2 kinase/fructose 2, 6 biphosphatase 3, PFKFB3 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	PFKFB3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	59,609
Uniprot ID	Q16875
Solubility (25°C)	6400ng/L
Sensitivity	10.25ng/L
Entrez	5209
Expiration Date	Please enquire.